

MICRO FILTER ASSEMBLY

4399

MICRO FILTER HOLDER ASSEMBLY, WITH GROUND GLASS JOINT

CAT. REF.	COMPONENTS	Qty Per Case
4399/47	Funnel for 47 size filter base with sintered disc for 34/35 size socket in the stem for 47 size filter spring clamp, 47 size 1000 ml. buchner flask with 34/35 cone	1



4396

MICRO FILTER HOLDER ASSEMBLY, WITH RUBBER CORK

CAT. REF.	COMPONENTS	QTY. PER CASE
4396/13	Funnel for 13 & 25 size filter base with sintered disc for 13 size filter spring clamp, 13 & 25 size 100 ml. buchner flask with cork	1
4396/25	Funnel for 13 & 25 size filter base with sintered disc for 25 size filter spring clamp, 13 & 25 size 100 ml. buchner flask with cork	1
4396/47	Funnel for 47 size filter base with sintered disc for 47 size filter spring clamp, 47 size 1000 ml. buchner flask with cork	1



4951

COMPLETE EXTRACTION APPARATUS

CAT. REF.	CAP. ML.	FLASH SIZE ML	QTY. PER CASE
4951/29	2000	5000	1
4951/32	5000	10000	1
4951/35	10000	20000	1

